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Please adhere to the following guidelines to ensure timely and accurate reports issued by Medical Research Institute

1. Please send all specimens through the laboratory of the respective hospital (except under special circumstances). A separate register may be maintained for this purpose in the laboratory.
2. The staff of the hospital laboratory is kindly requested to ensure that,
 - a. Whether the test is done at MRI before sending the samples. Contact MRI over the telephone if needed.
 - b. The correct specimen type is collected for the investigation. If not, inform the ward and request for properly collected specimen. (Eg: Tissue, Blood, Urine etc)
 - c. The specimen is collected into a proper container. If not, inform the ward and request for properly collected specimen. (Eg: Acid wash bottles for serum Ca⁺)
 - d. Separate specimens and request memos are sent for investigation at different laboratories at MRI. Eg: Sending a single specimen for both bacteriology and virology sections should be avoided.
 - e. Request memos are filled properly (if not inform the ward & request to complete before sending).
 - f. The address of the hospital is written in all request memos. Don't write abbreviations like THK. Place the rubber stamp of the hospital
 - g. The specimens are properly labeled as below.
 - h. A serial number is put both in the specimen label and the request (this may prevent the undue delay in sorting out specimens at the MRI specimen counter).
 - i. All specimens are entered in a delivery book according to the serial numbers as in the specimen.
 - j. Take protective measure to prevent the leakage of specimen during transport.
3. All the reports will be sent to the head of the laboratory of respective hospitals where a register is maintained for the delivery of laboratory reports to the wards. (A report issuing counter next to the specimen counter has been kept open for this task at MRI)
4. Transport media and special container for specimen collection could only be collected during office hours. (Weekdays 9.00 a.m. to 3.00 p.m., Saturday 9.00 a.m. to 11.30 am.
5. For more details please contact respective laboratories of the MRI during working hours and Resident MLT after hours.

Telephone: 0112693532 – 4, 0112633527 Hot line: 0112633536

Specimen label

Blood for Dengue antibody	
Name	: L.H. Sandali
BHT No	: 685163
Wd No	: 16
Hospital	: BH, Negombo

N.B.

Please take measures to bring the above contents to the notice of all relevant unit staff.