

DEPARTMENT OF VIROLOGY MEDICAL RESEARCH INSTITUTE DR DANISTER DE SILVA MAWATHA, COLOMBO 08

Tel: 011 269 3532-4 (Ext. 452)

Fax: 011 269 1495

Hepatitis B and C Real-time PCR

Sample collection / storage / transport

- 5 ml of blood sample collected into a disposable EDTA tube should be sent to MRI in ice
 (4⁰C) as early as possible (maximum delay 48 hours)
- The sample should be stored at 4°C until transport
- Sample collection can be conducted at MRI for patients between 8 am 12 pm on week days (Monday to Friday)

Important points

- Hepatitis B/C real-time PCR is an expensive test and should be ordered only when there is a real indication.
- The test should always be authorized by a Consultant.
- Serology results should always be available before requesting this test (Special situations should be discussed with the Consultant Virologist).
- 2 ml of blood should be sent in a plain bottle for the confirmation of HCV Ab result along with the sample for HCV PCR.